

EXAMVISION™  
Magnification  
Ergonomics



The perfect loupe



Magnifications available 2.3x 2.8x 3.5x 4.2x



## No compromises

To create the perfect loupe, we have developed a unique multi-lens system. The Galilean through-the-lens (TTL) design brings the loupe closer to the eye, increasing both field of view and depth of field.

Our optics and lens coatings are manufactured to our specifications in Germany, to give optimal light and true colour transmission.

## Custom made for you

Our optical specialists – usually optometrists – visit you to demonstrate the loupes.

They can assess how you work, and measure relevant factors to determine your specific needs.

They will carefully adjust a set of demonstration loupes, so you can experience how our loupes and headlights improve your working day.

Your loupes will be individually aligned and adjusted at our factory in Denmark to meet your specifications.



## Working in comfort

We have created lightweight, pure titanium frames. They are perfectly balanced, to ensure the loupes sit comfortably.


Naturally, our loupes help to improve your ergonomic position, so reducing neck and back strain while you work. It is possible to achieve an even more upright posture with our optional Prisma Lens System™. Let our optical specialists advise you.

## Strengthen your position

You want to give your clients the best possible treatment. So do we.

We offer a range of Halogen, HID and LED headlights tailor-made for our loupes, and a selection of accessories that help you work in comfort and safety.





**[www.exam-vision.com](http://www.exam-vision.com)**

Each pair of ExamVision loupes is custom made at our factory on the small island of Samsø in Denmark.

Founded by Jørn Rønvig (Rønvig Dental Mfg) and optometrist Kim Jensen in 2001, ExamVision combines unique optical systems and the latest in technology to produce the perfect TTL loupe.

*Danish design means it doesn't hurt to look good in loupes!*

**T: +45 87 92 12 10**  
**E: [mail@exam-vision.com](mailto:mail@exam-vision.com)**

